



Deutsche
Forschungsgemeinschaft



JPND
research

EU Joint Programme – Neurodegenerative Disease Research




Organization of neurophysiological and motion data for open and reproducible science

4th december 2025

iMOVE - Lübeck




Julius Welzel






I-MOVE
2025




Carl von Ossietzky
Universität
Oldenburg



What we can learn from animal research




Position only Direction only Position + Direction





Behaviour
labels

Time
labels


Neural data
(N)



Neurobehavioral “phenotyping”



Levels of „mobile“ brain recordings



- Data collection
(Hands-on session: Martin & Melanie)
- EEGManySteps
(Poster 13: Klaus)
- Gait related EEG
(Poster: 38 Lila)




The Brain Imaging Data Structure (BIDS)

What is BIDS?



Dataset info

Motion



Aims of BIDS




The Structure (The Container)

- **Hierarchy:** Standardized Subject → [Session] → modality folders.
- **Predictable Naming:** Key-value pairs allow automated parsing.
- **Simple Formats:** Uses non-proprietary standards (.edf, .tsv).


Metadata (The Context)

- **JSON Sidecars:** Every data file is paired with a .json file containing critical metadata details (e.g. sampling rate).
- **Machine-Readable:** Allows software (fMRIPrep, fieldtrip) to easily read data.
- **Beyond “Brain” Imaging:** dataset, participant, events descriptions



The Past, Present, and Future of BIDS

BIDS history timeline



My favourite two BIDS principles



Requirement levels

Level of requirement **for** files and **within** files:

1. **REQUIRED**

2. **RECOMMENDED**

3. **OPTIONAL**

```
├── README.md
├── dataset_description.json
├── participants.tsv
├── sub-001
│   ├── ses-01
│   │   ├── eeg/
│   │   └── motion/
│   └── ses-02
└── sub-002
```

participant_id	age	group
sub-001	34	C
sub-002	12	P
sub-003	33	M


Directories and file naming structure

- Data for each subject are organized in subdirectories labeled "**sub-<label>**"
- Each subject contains subdirectories for different data (type) types
- For a data file that was collected from a given **subject**, the file name **MUST** begin with the string **sub-<label>** and **MUST** contain a task-<label>

```
├── README.md
├── dataset_description.json
├── participants.json
├── participants.tsv
└── sub-001
    ├── eeg/
    │   └── sub-001_task-GoNoGo.eeg
```

Dataset info

EEG



- README.md
- dataset_description.json
- participants.json
- participants.tsv
- └ sub-001
 - ses-01
 - eeg/
 - motion/
 - ses-02
- └ sub-002

Modality agnostic data

dataset_description.json

├── README.md

├── **dataset_description.json** →

├── participants.json

├── participants.tsv

└── sub-001

 ├── ses-01

 └── ses-02

└── sub-002

```
{  
  "Name": "The mother of all experiments",  
  "BIDSVersion": "1.4.0",  
  "DatasetType": "raw",  
  "License": "CC0",  
  "Authors": [ "Paul Broca", "Carl Wernicke" ],  
  "Funding": [ „NIH F37823MFH1" ],  
  "EthicsApprovals": [ „HRPO (Protocol AR0928" ],  
  "DatasetDOI": "10.0.2.3/dfjj.10"  
}
```


Participants.tsv/.json

- └─ README.md
- └─ dataset_description.json
- └─ **participants.tsv**
- └─ **participants.json**
- └─ sub-001
 - └─ ses-01
 - └─ ses-02
- └─ sub-002

participants.tsv

participant_id	age	group
sub-001	34	C
sub-002	12	P
sub-003	33	M

participants.json

```
"age": {  
  "Description": "age of the  
  participant", "Units": "years"  
},  
„group“: {  
  "Description": „assigned group“,  
  "Levels": {  
    „C“: „control“,  
    „M“: „medication“,  
    „P“: „placebo“,  
  }  
}
```

Dataset info

EEG




- |— README.md
- |— dataset_description.json
- |— participants.json
- |— participants.tsv
 - └ sub-001
 - |— ses-01
 - |— eeg/
 - |— motion/
 - └ ses-02
 - └ sub-002

Modality specific data

Overview | eeg

└─ sub-001
└─ ses-01
└─ eeg/



```
{  
  "TaskName": "TASKNAME",  
  "SamplingFrequency": 1000,  
  "SoftwareFilters": "n/a",  
  "EEGChannelCount": 4,  
  "EOGChannelCount": 1,  
  "EEGReference": "placed on Cz",  
  "PowerLineFrequency": 50  
}
```

└─ sub-001_ses-01_task-<label>_eeg.<extension>

└─ sub-001_ses-01_task-<label>_eeg.json

└─ sub-001_ses-01_task-<label>_channels.tsv

└─ sub-001_ses-01_task-<label>_electrodes.tsv

└─ sub-001_ses-01_task-<label>_electrodes.tsv

name	type	units	status	status_description
CP5	EEG	microV	good	n/a
FC5	EEG	microV	bad	high freq noise
FC1	EEG	microV	good	n/a
C3	EEG	microV	good	n/a
VEOG	EOG	microV	good	n/a

name	x	y	z	impedance
CP5	-0.77	-0.30	0.57	8
FC5	-0.77	0.30	0.57	12
FC1	-0.29	0.31	0.91	2
C3	-0.59	0.00	0.81	5
VEOG	n/a	n/a	n/a	n/a

```
{  
  "EEGCoordinateSystem": "T1w",  
  "EEGCoordinateUnits": "mm",  
  "AnatomicalLandmarkCoordinates": {  
    "LPA": [-0.067, 1.736e-09, -3.844e-09],  
    "NAS": [-4.11e-09, 0.091, -4.541e-10],  
    "RPA": [0.064, -6.435e-09, -4.566e-09]  
  },  
  "AnatomicalLandmarkCoordinateSystem": "T1w",  
  "AnatomicalLandmarkCoordinateUnits": "mm",  
  "IntendedFor": "sub-01_T1w.nii.gz"  
}
```

Metadata | eeg.json

REQUIRED

TaskName, SamplingFrequency,
PowerlineFrequency, EEGReference,
SoftwareFilters

RECOMMENDED

TaskDescription, Manufacturer,
EEGChannelCount,
EOGChannelCount, ...

OPTIONAL

-

```
{  
  TaskName: Seeing stuff,  
  SamplingFrequency: 2400,  
  Manufacturer: Brain Products,  
  ManufacturersModelName: BrainAmp DC,  
  EEGChannelCount: 87,  
  EOGChannelCount: 2,  
  TriggerChannelCount: 1,  
  PowerLineFrequency: 50,  
  EEGPlacementScheme: 10-20,  
  EEGReference: single electrode placed on FCz,  
  SoftwareFilters: {  
    Anti-aliasing filter: {  
      half-amplitude cutoff (Hz):  
        500,  
      Roll-off: 6dB/Octave  
    }  
  }  
}
```

Metadata | channels.tsv

REQUIRED

name, type, units

RECOMMENDED

-

OPTIONAL

description, sampling_frequency

name	type	units	description	reference	status	status_description
VEOG	VEOG	uV	left eye	VEOG-, VEOG+	good	n/a
FDI	EMG	uV	left first dorsal interosseous	FDI-, FDI+	good	n/a
Cz	EEG	uV	n/a	REF	bad	high frequency noise
UADC001	MISC	n/a	envelope of audio signal	n/a	good	n/a


Metadata | electrodes.tsv

name	x	y	z	type	material	color
VEOG+	n/a	n/a	n/a	cup	Ag/AgCl	blue
VEOG-	n/a	n/a	n/a	cup	Ag/AgCl	white
FDI+	n/a	n/a	n/a	cup	Ag/AgCl	red
FDI-	n/a	n/a	n/a	cup	Ag/AgCl	red
GND	-0.07	0.00	-0.070	clip-on	Ag/AgCl	pink
Cz	0.00	0.07	0.06	cup	Ag/AgCl	yellow
REF	-0.07	-0.02	-0.01	cup	Ag/AgCl	grey

Column name	Requirement Level	Data type	Description
name	REQUIRED	string	Name of the electrode contact point. Values in <code>name</code> MUST be unique. This column must appear first in the file.
x	REQUIRED	number	Recorded position along the x-axis. This column must appear second in the file.
y	REQUIRED	number	Recorded position along the y-axis. This column must appear third in the file.
z	REQUIRED	number or "n/a"	Recorded position along the z-axis. This column must appear fourth in the file.
type	RECOMMENDED	string	Type of the electrode (for example, cup, ring, clip-on, wire, needle). This column may appear anywhere in the file.
material	RECOMMENDED	string	Material of the electrode (for example, <code>Tin</code> , <code>Ag/AgCl</code> , <code>Gold</code>). This column may appear anywhere in the file.
impedance	RECOMMENDED	number	Impedance of the electrode, units MUST be in <code>kOhm</code> . This column may appear anywhere in the file.
Additional Columns	OPTIONAL	n/a	Additional columns are allowed if they are defined in the associated metadata file.

Dataset info

Motion



- └─ README.md
- └─ dataset_description.json
- └─ participants.json
- └─ participants.tsv
- └─ sub-001
 - └─ ses-01
 - └─ eeg/
 - └─ **motion/**
 - └─ ses-02
- └─ sub-002

Modality specific data

Files | motion

└─ sub-001/
 └─ ses-01/
 └─ eeg/
 ...
 └─ motion/

- └─ sub-001_ses-01_task-<label>_tracksys-<label>_motion.tsv
- └─ sub-001_ses-01_task-<label>_tracksys-<label>_motion.json
- └─ sub-001_ses-01_task-<label>_channels.tsv

0,2634511	0,092295	0,0086682	0,9305117	0,690106	0,8098815
0,694520	0,1918243	0,8437273	0,3975710	0,885496	0,8952724
0,0766395	0,2587211	0,5434792	0,2822837	0,2789791	0,2326254
0,5779993	0,0456141	0,0490745	0,9408899	0,1533421	0,6683652
0,054556	0,7915927	0,5871733	0,4669577	0,9754468	0,0480541
0,966024	0,1962834	0,7114406	0,3389448	0,7194495	0,4384892
0,984172	0,5079461	0,1180168	0,7966978	0,1753768	0,4886533
0,9883907	0,1557346	0,8002013	0,6334882	0,7526906	0,8529441
...

```
{  
  "SamplingFrequency": 60,  
  "SamplingFrequencyEffective": 60.19,  
  "TaskName": "BIDS Motion fictive example",  
  "TrackingSystemName": "imu1",  
  "TaskDescription": "walking and talking",  
  "MotionChannelCount": 6,  
  "SubjectArtefactDescription": "n/a",  
  "TrackedPointsCount": 2,  
  "ACCELChannelCount": 3,  
  "GYROChannelCount": 3,  
  "Manufacturer": "BWSensing",  
  "ManufacturersModelName": "BW-imu600",  
}
```

name	component	type	tracked_point	units	placement
imu1_rf_acc_x	x	ACCEL	rf	m/s ²	right_foot
imu1_rf_acc_y	y	ACCEL	rf	m/s ²	right_foot
imu1_rf_acc_z	z	ACCEL	rf	m/s ²	right_foot
imu1_rf_gyro_x	x	GYRO	rf	rad/s	right_foot
imu1_rf_gyro_y	y	GYRO	rf	rad/s	right_foot
imu1_rf_gyro_z	z	GYRO	rf	rad/s	right_foot

Metadata | motion.json

REQUIRED

TaskName, SamplingFrequency

RECOMMENDED

TaskDescription, RecordingDuration,
MotionChannelCount, <type>ChannelCount,
SoftwareFilters, ...

OPTIONAL

Manufacturer, RecordingSoftware, ...

```
{  
  "TaskName": "BIDS Motion fictive example",  
  "SamplingFrequency": 60,  
  "TrackingSystemName": "IMU Right Hand",  
  "TaskDescription": "walking and talking",  
  "InstitutionAddress": "Fictive address",  
  "InstitutionName": "Fictive Institution",  
  "MotionChannelCount": 18,  
  "RecordingDuration": 4667.641106,  
  "SubjectArtefactDescription": "n/a",  
  "TrackedPointsCount": 2,  
  "ACCELChannelCount": 6,  
  "GYROChannelCount": 6,  
  "MAGNChannelCount": 6,  
  "Manufacturer": "BWSensing",  
  "ManufacturersModelName": "BW-IMU600",  
}
```

Metadata | channels.tsv

REQUIRED

name, component, type, tracked_point, units

RECOMMENDED

placement, reference_frame

OPTIONAL

status, status_description, sampling_frequency


name	component	type	tracked_point	units	reference_frame
t1_acc_x	x	ACCEL	LeftFoot	m/s^2	global
t1_acc_y	y	ACCEL	LeftFoot	m/s^2	global
t1_acc_z	z	ACCEL	LeftFoot	m/s^2	global
t1_gyro_x	x	GYRO	LeftFoot	rad/s	global
t1_gyro_y	y	GYRO	LeftFoot	rad/s	global
t1_gyro_z	z	GYRO	LeftFoot	rad/s	global
...					

Restricted keyword list for channel type

Restricted keyword list for column `type` in alphabetic order. No

Keyword	Description
ACCEL	Accelerometer channel, one channel for each spatial axis. <code>*_channels.tsv</code> file (x, y, or z).
ANGACCEL	Angular acceleration channel, one channel for each spatial axis. <code>*_channels.tsv</code> file (x, y, or z).
GYRO	Gyrometer channel, one channel for each spatial axis. <code>*_channels.tsv</code> file (x, y, or z).
JNTANG	Joint angle channel between two fixed axis belonging to two proximal and distal bodypart in deg.

Overview | motion



BIDS Starterpack

Website

[<https://bids-specification.readthedocs.io/en/stable/>]

Tutorials

[MNE, fieldtrip, EEGLab]

Validator

[<https://bids-standard.github.io/bids-validator/>]

Community

[https://juliuswelzel.github.io/igor_bids_consultation/]



MNE-BIDS 0.16.0 document...[News](#)[Install](#)[Use](#)[API](#)[CLI](#)[Contribute](#)

BIDS V

IGOR BIDS consultations hourMeetingsMaterials

IGOR BIDS consultations hour

It's a match!



What is the IGOR BIDS Consultation Hour?

We're are happy to offer the IGOR BIDS Consultation Hour, a monthly forum designed to help researchers like you navigate the Brain Imaging Data Structure (BIDS) format. During these sessions, you can connect with experienced BIDS users to get individualized feedback and solutions for your specific BIDS data conversion challenges.

What can I expect?

- **Personalized guidance:** Our experienced researchers, [Julius](#) (EEG & motion data expertise) and [Lennart](#) (MRI expertise), will be available to answer your questions and provide tailored advice on converting your datasets to the BIDS format.
- **Open forum:** Feel free to bring any BIDS-related questions you have, big or small. We're here to help you overcome conversion hurdles and ensure your data adheres to BIDS best practices.
- **Recurring sessions:** The IGOR BIDS Consultation Hour takes place on the 4th Friday of each month at 12:00 PM (noon) CET.

read BIDS datasets

Convert MNE sample data to BIDS format


interactive data inspection and bad channel selection

Convert EEG data to BIDS format


On this page
[It's a match!](#)
What is the IGOR BIDS Consultation Hour?
What can I expect?
How can I participate?
Upcoming Session:

BIDS in the bigger picture

- Building a sustainable research environment
- Connect rich context information with neural recordings



Conclusions



Thanks



Sein Jeung & Team






Lara Papin

Hands-On session

https://juliuswelzel.github.io/motion_bids_materials_jw



How to combine EEG and motion data?

